



University of Colorado
Anschutz Medical Campus

Comparison of Clinic-Based Biomechanical Measures During Walking to Laboratory Measures after Total Knee Arthroplasty

Bade M^{1,2}, Edwards E¹, Juan J¹, Sapinsley Z¹, Thomsen PB¹, Peters A^{1,2}, Zeni J Jr⁴, Cheuy VA⁵, Christiansen CL^{1,2}, Stevens-Lapsley J^{1,2}

¹University of Colorado Anschutz Medical Campus, Aurora, CO

²VA Eastern Colorado Geriatric Research, Education, and Clinical Center (GRECC), VA Eastern Colorado Health Care System, Aurora, CO

³Rutgers University, school of Health Professions, Newark, NJ

⁴VA Rocky Mountain Mental Illness, Research, Education and Clinical Center (MIRECC), VA Eastern Colorado Healthcare System, Aurora, CO

Scan me to learn more





PURPOSE

Clinically feasible methods to assess movement biomechanics post total knee arthroplasty are needed. However, the relationship between measurements collected clinically to those collected in a laboratory setting has not been examined. Therefore, the purpose of this study was to compare clinic-based measures of gait biomechanics using insole sensors collected during a physical therapy session to laboratory-based motion capture measurements 10-weeks after TKA.

ACKNOWLEDGEMENTS

This research was funded : NIH R01 AG056585, NIH/NCATS UL1-TR001082, NIH F32-AG066274, NCATS KL2 TR001870, Eastern Colorado Geriatric Research Education and Clinical Center, Academy of Orthopaedic Physical Therapy Career Development Award

Insole pressure measurement systems may have **clinical utility** as a **less-costly alternative** to **motion capture systems** for investigating between-limb kinetic symmetry.



Figures 1 and 2. Laboratory-Based Analysis Using Motion Capture and Embedded Force Plates of Self-Selected Walking Speed at 10 weeks after TKA



Figure 3. Clinic-Based Analysis Using Single Sensor Insoles of Self-Selected Walking Speed at 10 weeks after TKA (image credit www.NovelUSA.com/loadsol)

	Gait Mean (SD)	Clinic Mean (SD)	Paire t-test p-value	Correlation (95% CI)	ICC (95% CI)	Bias (95% LoA)
Involved Loading Rate (BW/s)	8.43 (3.7)	6.781 (2.1)	0.0006	0.34* (0.10, 0.54)	0.451 (0.09, 0.67)	1.65 (-5.32, 8.63)
Uninvolved Loading Rate (BW/s)	10.31 (3.9)	7.499 (2.4)	<0.0001	0.28** (0.03, 0.50)	0.405 (0.008, 0.64)	2.81 (-4.97, 10.6)
Loading Rate Ratio	0.848 (2.7)	0.927 (0.19)	0.07	-0.06 (-0.30, 0.20)	-0.115 (-0.86, 0.33)	-0.079 (-0.74, 0.58)
Involved GRF (N/BW)	1.05 (0.09)	0.980 (0.08)	<0.0001	0.22 (-0.04, 0.44)	0.353 (-0.08, 0.61)	0.066 (-0.14, 0.27)
Uninvolved GRF (N/BW)	1.09 (0.09)	1.06 (0.22)	0.34	0.04 (-0.21, 0.29)	0.059 (-0.57, 0.44)	0.028 (-0.42, 0.48)
GRF Ratio	0.962 (0.04)	0.940 (0.10)	0.14	-0.01 (-0.26, 0.24)	-0.016 (-0.69, 0.39)	0.022 (-0.20, 0.25)
Involved Impulse (BW*s)	0.538 (0.04)	0.609 (0.17)	0.002	0.14 (-0.12, 0.38)	0.120 (-0.47, 0.47)	-0.070 (-0.41, 0.27)
Uninvolved Impulse (BW*s)	0.567 (0.05)	0.617 (0.08)	<0.0001	0.27** (0.02, 0.49)	0.389 (-0.02, 0.63)	-0.050 (-0.22, 0.12)
Impulse Ratio	0.952 (0.06)	0.980 (0.14)	0.20	-0.28** (-0.50, -0.03)	-0.517 (-1.53, 0.09)	-0.028 (-0.36, 0.30)

Table 1. Summary Statistics, Correlations, ICC (3,k), Bias (n=61)
*p<0.01; **p<0.05. Abbreviations: BW, bodyweight




Figure 4. Bland-Altman Plot of Involved Loading Rate




Figure 5. Bland-Altman Plot of Ground Reaction Force Ratio

PARTICIPANTS

61 participants (aged 65±8 years)

METHODS

Subjects underwent biomechanical assessments while walking 10-weeks after TKA: 1) using single sensor insoles during an intervention session (clinic-based assessment) and 2) using an 8-camera motion capture system and embedded force plates (laboratory-based assessment). Average vGRF, impulse, and loading rate for the surgical limb, uninvolved limb, and between-limb symmetry ratio (surgical/uninvolved) were collected.




Figure 6: Vertical Ground Reaction Force Curve Over 0-100% of the Gait Cycle

RESULTS

Symmetry ratios and uninvolved-limb vGRF were not different between clinic and laboratory-based assessments (p>0.05) but were different for surgical limb vGRF, impulse and loading rate (p<0.05) and uninvolved limb impulse and loading rate (p<0.05).

CLINICAL RELEVANCE

Clinic-based symmetry ratios demonstrated greater accuracy than individual limb measurements of vGRF, impulse, and loading rate with laboratory-based measurements of walking 10-weeks after TKA. This may be due to ratios having the advantage of correcting for systematic differences in magnitude (e.g. insoles forces being lower on both surgical and uninvolved limbs).